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Edition 02.15

**Laboratory Control Unit LCU,  
Laboratory Safety Valve VCL,  
Central Shut-off Device VCC**

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**⚠️ WARNING**

**Explosion hazard if gas installation tasks  
are carried out improperly!**

Incorrect installation, adjustment, modification,  
operation or maintenance may cause injury or material  
damage.

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**Read the instructions before use.**

This unit must be installed in accordance with the  
regulations in force.

Note DVGW Code of Practice G 621.

All the work set out in these operating instructions may  
only be completed by trained and authorized personnel!

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**Contents**


1	Declaration of Conformity .....	2
2	Product Liability and Warranty .....	2
2.1	Scope of Delivery .....	2
2.2	Intended Use .....	2
2.3	Warranty .....	2
2.4	User's Inspection .....	2
3	Installation .....	3
3.1	Requirements .....	3
3.1.1	System Volume .....	3
3.2	Installation of VCL/VCC .....	3
3.3	Installation of LCU .....	4
3.3.1	Plastic housing LCU..A .....	4
3.3.2	Plastic Housing of LCU..E .....	5
3.3.3	Plastic Housing for Concealed Installation of LCU..U .....	5
3.4	Wiring .....	5
3.4.1	Preparing the Wiring .....	5
3.4.2	VCL (V1, V2), VCC (V1, V2) .....	6
3.4.3	Venting valve VBY (VE) of VCL .....	6
3.4.4	Pressure sensor (PIC) of VCL .....	6
3.4.5	Electrical connection: LCU with Key Operated Switch .....	7
3.4.6	LCU without Key-Operated Switch, with Daughterboard .....	8
3.5	Assembly After Wiring .....	8
4	Tightness Test .....	9
5	Commissioning .....	9
5.1	Preparatory Tasks for Commissioning .....	9
5.2	Performing Initial Commissioning .....	10
6	Operation .....	11
6.1	Open Central Shutoff Valve .....	11
6.2	Consumer Startup .....	12
6.3	Power Failure .....	12
6.4	Timeout .....	12
6.5	System Shutdown .....	12
7	Parameter Indication .....	13
7.1	Factory Settings .....	13
8	Maintenance .....	14
9	Technical Advice .....	14
10	Accessories .....	15
11	Technical Data .....	16

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
## 1 Declaration of Conformity

### LCU

We, the manufacturer, hereby declare that the product LCU, marked with product ID number CE-0085BS5046, complies with the essential requirements of the following Directives:

- 90/396/EEC in conjunction with EN 13611, EN 13611-A1, EN 1643, EN 298,
- 2006/95/EC in conjunction with the relevant standards,
- 2004/108/EG in conjunction with the relevant standards,
- 2009/142/EG in conjunction with the relevant standards.
- The correspondingly marked products comply with the tested type sample. Comprehensive quality assurance is guaranteed.

### VCL/VCC

The laboratory safety valve VCL and the central shut-off device VCC is a combination of two VAS gas solenoid valves. We, the manufacturer, hereby declare that the product VAS, marked with product ID No. CE-0063BO1580, complies with the essential requirements of the following Directives:

- 90/396/EEC in conjunction with EN 161, EN 88, EN 126, EN 12067-1, EN 1854,
- 2006/95/EC in conjunction with the relevant standards,
- 2004/108/EG in conjunction with the relevant standards,
- 2009/142/EG in conjunction with the relevant standards.

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
## 2 Product Liability and Warranty

This product carries the CE-mark, which means that it satisfies the essential requirements laid down in the EC Directives.

We are the manufacturer of this product:

**Gastechnik Kirchner GmbH & Co. KG**

- Gröninger Weg 7, D-74379 Ingersheim
- Phone +49 (0)7142 9191-30
- Fax +49 (0)7142 9191-40
- info@gastechnik-kirchner.de
- www.gastechnik-kirchner.de

### 2.1 Scope of Delivery

1. Laboratory control unit LCU
2. 2 x keys for key-operated switch, not applicable to LCU..M / LCU..S
3. 5 x M16 cable gland
4. 1 x M20 cable gland
5. VCL comprising:
  - double solenoid valve VCS
  - venting valve VBY
  - pressure sensor with plug and 5 m connection cable
6. VCC comprising:
  - double solenoid valve VCS

### 2.2 Intended Use LCU

Laboratory control unit LCU for fully automatic checking of the downstream installation for leaks, tight closing of valves and adequate inlet pressure before every system start-up.

### VCL

Laboratory safety valve VCL for safeguarding gas or air on various appliances.

### VCC

Central shut-off device VCC as shut-off valve upstream of VCL valve.

### 2.3 Warranty

Our Standard Terms and Conditions of Sale effective at the time shall apply. Agreements diverging from these Standard Terms and Conditions do not restrict the legal rights of the buyer.

Any warranty exceeding the above provisions shall require a contractual form and shall exclude component-related vandalism, software updates and consumables.

### 2.4 User's Inspection

Immediately upon receipt, the goods must be checked for completeness and potential damage in transit. Notice of any such damage must be given immediately.

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

- ☐ Natural gas
- ☐ LPG (gaseous)
- ☐ Air

#### ⚠ CAUTION

The gas must be dry in all temperature conditions and must not contain condensate.

For mains voltage, electrical power consumption, ambient temperature, type of protection, inlet pressure  $p_e$  and installation position, see type label.

#### ⚠ CAUTION

Do not exceed the maximum inlet pressure  $p_e$  of 100 mbar. Non-compliance can damage the pressure sensor or the valve.

##### LCU

Only to be used in grounded mains.

##### VCL/VCC

##### Opening times:

Fast opening:  $\leq 0,5$  s

##### Closing time:

Fast closing:  $< 1$  s

The inlet pressure  $p_e$  can be measured using a pressure test fitting, see section "Accessories", page 15.

#### ⚠ CAUTION

The solenoid actuator heats up during operation. Surface temperature  $< 85^\circ\text{C}$  ( $185^\circ\text{F}$ ) acc. to EN 60730-1 at an overvoltage of 10% and an ambient temperature of  $25^\circ\text{C}$  ( $77^\circ\text{F}$ ). Use a heat-resistant cable rated at  $85^\circ\text{C}$  ( $185^\circ\text{F}$ ).

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### 3.1.1 System Volume

The system volume is limited to maximum 60 l ( $= \wedge 105$  m piping with nominal diameter DN 25). The minimum system volume downstream of the VCL is 0.2 l ( $= \wedge 0.20$  m piping with nominal diameter DN 25).

System volume VP is calculated from the valve volume VV, added to the volume of the piping VR for each additional meter in length L ( $VP = VV + L \times VR$ ). If the piping volume (on-site installation, system expansion, conversion), the inlet pressure range or the gas type changes, or if the laboratory control unit LCU or the laboratory safety valve VCL are replaced, commissioning must be carried out again by the fitter/authorized personnel (see section "Commissioning", page 9).

Ventile Valves	Ventilvolumen VV [ l ] Valve volume VV [ l ]	DN	Volumen der Rohrleitung pro Meter VR [ l/m ] Volume of piping per meter VR [ l/m ]
VCL 115 / VCC 115	0,05 / 0,05	15	0,2 / 0,2
VCL 120 / VCC 120	0,05 / 0,05	20	0,3 / 0,3
VCL 125 / VCC 125	0,05 / 0,05	25	0,5 / 0,5

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### 3.2 Installation of VCL/VCC

#### Fitting position of VCL/VCC

Black solenoid actuator in vertical or oblique to horizontal position, but not upside down.

Make sure that the distance from potential obstacles is always at least 20 mm in any direction.

Sealing material and thread cuttings must not be allowed to get into the valve housing and the sealing surfaces.

Ensure that there is sufficient space for installation and adjustment.

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
#### NOTICE

To keep the laboratory safety system free of dirt, we recommend that a filter ( $50 \mu\text{m}$ ) be installed upstream of the system.

A soiled laboratory safety system may lead to incorrect measurements or cause complete system failure.

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
**⚠ CAUTION**

Do not clamp the unit in a vice.  
Keep it in position with a suitable open-end wrench applied to the octagon of the flange.  
Non-compliance may cause housing damage: risk of external leakage.

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
**⚠ CAUTION**

Do not hold device by the solenoids!  
Grasp valve body!

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**3.3 Installation of LCU**

Fitting position of LCU: horizontal


- ☞ Mount laboratory control unit LCU outside of or inside the laboratory/training room.
- ☞ Choose a location where the laboratory control unit LCU will be kept free from vibrations.
- ☞ Housing must be fitted level (avoid mechanical stress).

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
**NOTICE**


For the cable leadthroughs: Use M20 and M16 cable glands delivered with the unit to obtain IP 54 type of protection.

Retaining screws  $\varnothing 3.5 \times 25$  mm are recommended to mount the bottom part of the housing.

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
### 3.3.2 Plastic Housing of LCU..E

Remove cover before wiring!  
Cut out the housing contour using the template (4, enclosed).  
Perform steps 1 to 7 as illustrated.  
Milled hole: d = 68 mm

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### 3.3.3 Plastic Housing for Concealed Installation of LCU..U

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### 3.4 Wiring Cable selection

#### LCU

Use mains cable and fuses suitable for the type of operation and complying with local regulations.  
3 x max. 1,5 mm<sup>2</sup>

#### VCL/VCC

Use a heat-resistant cable rated at 85°C (185°F).

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V1 and V2:  
4 x max. 1.5 mm<sup>2</sup>

Venting valve VBY (VE):  
3 x max. 1.5 mm<sup>2</sup>

Pressure sensor (PIC):  
LIYCY, max. 2 x max. 0.75 mm<sup>2</sup>, shielded and shield connected to the LCU (not to the plug).

#### If fitted:

LCU to VCC (individual valve EV):  
3 x max. 1.5 mm<sup>2</sup>


LCU to Emergency Stop button GNA:  
3 x max. 1.5 mm<sup>2</sup>

LCU to building services management system:  
Operating signal:  
2 x max. 1.5 mm<sup>2</sup>

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### 3.4.1 Preparing the Wiring

- ☞ Disconnect the system from the power supply.
- ☞ Shut off the gas supply.
- ☞ Wiring acc. to EN 60204-1.

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
- 1 = N
- 2 = L<sub>V1</sub>
- 3 = L<sub>V2</sub>


VCC connections

- 1 = N
- 2 = L<sub>V1/V2</sub> (internally bridged with jumper)
- 3 = not assigned

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


#### ⚠ CAUTION

Do not rotate the socket attached to the pressure sensor (see red circle mark)! Non-compliance can damage the sensor's internal wiring.

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
Do not shorten the connection cable. Just coil up any excess cable length, since the sensor cable is supplied preassembled and tested.

If the socket is to be redirected, the socket insert must be rotated accordingly.

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
**CAUTION**

A short-circuit on the valve connections or on the pressure sensor damages internal fine-wire fuses. These fuses cannot be replaced. External voltage applied to terminals 1 and 2 will destroy the LCU.


Wire the LCU as shown in the connection diagram.

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**3.4.5 Electrical connection:  
LCU with Key Operated Switch**


- BB** Anzeige  
Display
- V1** Eingangsseitiges Gasventil V1 im Labor-Sicherheitsventil VCL  
Inlet gas valve V1 in the laboratory safety valve VCL
- V2** Ausgangsseitiges Gasventil V2 im Labor-Sicherheitsventil VCL  
Outlet gas valve V2 in the laboratory safety valve VCL
- VE** Entlüftungsventil VBY  
Venting valve VBY
- EV** Gasventil V1/V2 im Einzelventil VCC  
Gas valve V1/V2 in individual valve VCC
- B** Konfigurierbares, potentialfreies Melderelais  
Configurable, potential-free transmitting relay
- S** Schlüsselschalter  
Key-operated switch
- NTA** Drucksensor  
Pressure sensor
- NTA** Not-Aus-Taster: **potentialfreier Öffner!**  
Emergency Stop button: **potential-free break contact!**
- BR** V1/V2 intern gebrückt in der VCC  
V1/V2 internally bridged in VCC

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### 3.4.6 LCU without Key-Operated Switch, with Daughterboard

- BB** Anzeige  
Display
- V1** Eingangsseitiges Gasventil V1 im Labor-Sicherheitsventil VCL  
Inlet gas valve V1 in the laboratory safety valve VCL
- V2** Ausgangsseitiges Gasventil V2 im Labor-Sicherheitsventil VCL  
Outlet gas valve V2 in the laboratory safety valve VCL
- VE** Entlüftungsventil VBY  
Venting valve VBY
- EV** Gasventil V1/V2 im Einzelventil VCC  
Gas valve V1/V2 in individual valve VCC
- B** Konfigurierbares, potentialfreies Melderelais  
Configurable, potential-free transmitting relay
- NTA** Not-Aus-Taster: **potentialfreier Öffner!**  
Emergency Stop button: **potential-free break contact!**
- BR** V1/V2 intern gebrückt in der VCC  
V1/V2 internally bridged in VCC

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
#### LCU..A

To fasten the cover, tighten the four screws with a Phillips screwdriver.

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


#### 4 Tightness Test

##### ⚠ CAUTION

Do not exceed the maximum inlet pressure  $p_e$  of 100 mbar. Non-compliance can damage the pressure sensor.

To be able to check the tightness in accordance with DVGW Code of Practice G 600, we recommend supplying nitrogen via a separate pressure test point, see section „Accessories“, page 15.  
Max. allowable pressure:  
 $p_{e \max} = 100 \text{ mbar}$

##### Requirements

- ☞ Disconnect the system from the power supply.
- ☞ Shut off the gas supply.

##### Connections for measuring sockets

(1)	(2) + (3)	(4)
Admission pressure side	Medium pressure side	Back pressure side

- ☞ In case of piping leaking: Check seal on valve or pipe system!
- ☞ VCL leaking: dismantle the VCL and return it to the manufacturer

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##### NOTICE

Back pressure side piping:  
Release test pressure before commissioning.

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#### 5 Commissioning

##### ⚠ CAUTION


Commissioning may only be carried out by authorized fitters/personnel. The unit may only be released for unsupervised operation once commissioning has been completed correctly.

##### NOTICE

Before commissioning, the entire upstream and downstream system must be vented and tight in order to ensure fault-free operation. Inadequately vented systems/piping cause malfunctions or availability problems in the VCC/VCL units.

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##### 5.1 Preparatory Tasks for Commissioning


- ☞ All consumers, e.g. side tables, must be connected and the sampling points fitted to them must be closed.
- ☞ Open ball valve, or upstream valve (VCC), upstream of the VCL.
- ☞ Vent the entire upstream piping.

- ☞ Switch on the LCU using the key-operated switch.
- ☞ „Version display“ appears briefly in the display.
- ☞ The message „Consumer off/Individual valve off“ is displayed.

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
##### Call up factory settings

With the Enter button (1) held down, operate first the „Consumer“ button (2) and then the „Individual valve“ button (3).

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- ⇒ The message “Factory settings” is displayed.
- ☞ Select gas type with corresponding inlet pressure range:
- ☞ Go to “Type of gas/Inlet pressure range” using the arrow button. The factory setting is 10–40 mbar, natural gas.
- ☞ Press the Enter button. The cursor at the end of the line starts flashing.
- ☞ Select the required gas type and inlet pressure range using the arrow button.
- ☞ Confirm by pressing the Enter button.

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
### NOTICE

In the event of deviating inlet pressures (even for brief pressure deviations), the settings for the on-site pressure regulator must be corrected, or alternatively an additional pressure regulator must be installed directly upstream of the laboratory safety valve VCL.

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#### Quit factory settings

- ☞ Go to the menu item “Factory default settings” using the arrow button.
- ☞ Confirm by pressing the Enter button.
- ☞ The following messages are displayed one after the other:
- ⇒ LCU self-test
- ⇒ Consumer off Individual valve off

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## 5.2 Performing Initial Commissioning

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
### NOTICE

All tapping points downstream of the VCL must be closed.

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

With the Enter button (1) held down, operate first the “Consumer” button (2) and then the “Individual valve” button (3).

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- ⇒ The display changes to “Sensor Test”.

- ☞ The Commissioning process starts automatically after the sensor test. The VCL opens and the downstream pipe system is filled with gas. The testing time is determined during commissioning.

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### NOTICE

With the sensor test, first all three valves open briefly, then the venting valve VBY opens and the blue LED in the plug lights up. After this, the blue LED on the first valve of the VCL and the blue LED in the plug of the venting valve VBY must flash 5 times alternately. Both LEDs on the VCL will then light up briefly. Afterwards, the blue LED on the second valve of the VCL and on the VBY light up until the testing time has been determined.

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- ⇒ The display changes to "Commissioning OK".
- ☞ Acknowledge commissioning by pressing the Enter button.
- ⇒ The display changes to "Consumer off/Individual valve off".
- ⇒ Commissioning is completed.
- ⇒ The system is ready for operation.
- ☞ Switch off the LCU using the key-operated switch. Remove the key and store in a safe place.
- ☞ Close the ball valve or the valve upstream of the VCL and disconnect the LCU from the power supply.

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#### NOTICE

If the piping volume (on-site installation), the inlet pressure range or the gas type changes, or if the laboratory control unit LCU or the laboratory safety valve VCL are replaced, commissioning must be carried out again by the fitter/authorized personnel.

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#### Commissioning protocol

The commissioning data can be entered by the fitter/authorized personnel here.

Inbetriebnahme vom Commissioning date		LCU Seriennummer LCU serial number	
Gasart Type of gas	Raum Room	Eingangsdruck p <sub>e</sub> (mbar) Inlet pressure p <sub>e</sub> (mbar)	
Ermittelter Betriebsdruck (mbar) Operating pressure (mbar) determined		Ermittelte Prüfzeit t <sub>E</sub> (ms) Test time t <sub>E</sub> (ms) determined	

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
## 6 Operation

### 6.1 Open Central Shutoff Valve

Open ball valve upstream of the central shut-off valve VCC.

#### CAUTION

All tapping points on the teacher's workbench must be closed.

- ☞ Ensure voltage supply to LCU.
- ☞ Switch on the LCU using the key-operated switch.
- ⇒ The following messages are displayed one after the other:  
Version, Consumer off/  
Individual valve off.
- ☞ Press button "Individual valve".
- ⇒ The message "Consumer off/  
Individual valve on"  
is displayed.
- ☞ The central shut-off valve VCC opens.
- ☞ Open the Bunsen burner and tapping point downstream of the central shut-off valve VCC.
- ☞ Ignite the Bunsen burner.
- ⇒ The teacher's workbench is ready for operation.

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## 6.2 Consumer Startup

- ☞ Open the ball valve upstream of the laboratory safety valve VCL (consumer).

### CAUTION

All tapping points on the consumer workbenches must be closed.

- ☞ Ensure voltage supply to LCU.
- ☞ Switch on the LCU using the key-operated switch.
- ⇒ *The display shows: Version, Consumer Off, Indiv. Valve Off.*

### NOTICE

If the central shutoff valve is located upstream, it must first be opened. See section "Factory Settings", page 13, installation according to G 621.

- ☞ Press the "Consumer" button.
- ⇒ *First all three valves open briefly, then the blue LED on the first valve of the VCL lights up alternately with the blue LED in the plug of the venting valve VBY, max. 4 times each. Then both LEDs on the VCL light up briefly.*

The following messages are displayed one after the other:

- Sensor test
- Tightness test
- Tightness test OK
- ☞ The VCL opens and releases the gas supply.
- ⇒ *The consumer workbenches are ready for operation.*
- ☞ Open the tapping points and Bunsen burners downstream of the laboratory safety valve VCL (consumer).
- ☞ Ignite the Bunsen burners.

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
## 6.3 Power Failure

After a power failure during operation, all valves close in less than 1 s.

Once the power supply has been switched on again, the LCU shows the default setting "Consumer off / Individual valve off".

## 6.4 Timeout

If the "Consumer" button is not pressed after the power supply and the laboratory control unit LCU have been switched on, the LCU automatically shows the program status "Timeout" after 24 h.

In order to be able to use the LCU again after this, the LCU must be switched off and on again using the key-operated switch.


The default setting "Consumer off / Individual valve off" is displayed.

## 6.5 System Shutdown

Upon operating the "Consumer" button or "Individual valve" button, all gas valves will close in less than 1 s.

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## 7.1 Factory Settings

The following factory settings can be called up, and  
partly changed, by the fitter/authorized personnel.

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### Description

- Countdown has been factory-  
preset to 0 h: Consumer valves  
remain open.
- Countdown can be changed in  
the “Factory default settings”  
menu to max. 10 h.  
10 min before the end of the  
countdown is reached, a  
warning signal will be sounded.  
The warning tone can be  
switched off by pressing the  
“Consumer” button and the  
countdown will start again as a  
result. If the tone is not  
switched off, the message  
“Please restart LCU” will be  
displayed.

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INSTALLATION  
NACH G621: 0  
  
INSTALLATION  
AS PER G621: 0

- Installation according to G 621:  
0: No query (default)  
1: If teaching valve is not open  
at Start Consumer,  
message “Error Open  
Teaching valve”,  
VCL does not open.  
2. Teaching valve opens  
automatically at Start  
Consumer.

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FUNKTIONSART  
MELDERELAIS: 0  
  
FUNCTIONAL MODE  
TRANSMITTING RELAY: 0

- Function of signaling relay  
terminals 7/8/9  
0: Operating message:  
VCL open (default)  
1: Operating message:  
LCU ready  
2. Operating message:  
LCU collective fault  
3. Message: Test in progress  
4. Message: Insufficient gas  
5. Message: EV not open  
(acc. to G 621). Only active  
if installation acc. to G 621,  
option 1.
- Gas type with inlet pressure  
range

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BETRIEBSDRUCK  
20MBAR  
  
OPERATING PRESSURE  
20 MBAR

- Operating pressure determined  
Pressure range 10–40 mbar:  
Default operating pressure  
20 mbar  
Pressure range 25–70 mbar:  
Default operating pressure  
50 mbar

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## NOTICE

If the current operating pressure lies outside the  
determined operating pressure (tolerance +/- 5 mbar),  
the sensor test will fail.

☞ Perform commissioning

☞ If incoming pressure  $p_e$ , varies strongly, a pressure  
regulator may need to be incorporated upstream.

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2.720 ZAEHLER  
LCU  
  
2720 LCU  
COUNTER

- Number of power-on operations  
for the laboratory control unit  
LCU

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
## 8 Maintenance

### NOTICE

To ensure smooth operation, check VCL/VCC for tightness and proper functioning once every year.

- ☞ Disconnect the system from the power supply.
- ☞ Shut off the gas supply.

### Cleaning the strainer

If the flow rate is correct, see section "Tightness Test", page 9.

If the flow rate has dropped, clean the strainer.

The valves must only be removed from the piping and re-installed on the inlet and outlet flange.

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## 9 Technical Advice

If malfunction occurs, the laboratory control unit LCU closes the gas valves. At the same time, a corresponding error message is displayed and a fault signal is emitted.

### Error 01/Error Sensor Test

Pressure sensor check has failed.

- ☞ The current operating pressure lies outside the determined operating pressure (tolerance +/- 5 mbar)
- ⇒ Perform commissioning
- ⇒ If incoming pressure  $p_e$  varies strongly, a pressure regulator may need to be incorporated upstream.
- ☞ The nozzle on venting valve VBY is clogged.
- ⇒ Remove the unit and return it to the manufacturer for inspection.
- ☞ The cables for the valves VCL and VBY have been reversed.
- ⇒ Check wiring plan and correct, see section "Wiring", page 5.

### NOTICE

During commissioning, first all three valves open briefly, then the venting valve VBY opens. The blue LED in the plug lights up. After this, the blue LED on the first valve of the VCL and the blue LED in the plug of the venting valve VBY must flash five times alternately. Both LEDs on the VCL will then light up briefly. Afterwards, the blue LEDs on the second valve of the VCL and on the VBY light up until the testing time has been determined.

- ☞ Faulty pressure sensor signal.
- ⇒ Open laboratory control unit LCU. The measuring instrument is to be connected in series between the „+“ or „-“ terminal and the disconnected cable.
- ⇒ Measure the current signal: correct current signal range 4–20 mA (4 mA = 0 mbar)
- ⇒ Check the pipeline.
- ⇒ If the measures described above do not help, remove the unit and return it to the manufacturer for inspection.

### Error Commissioning

Commissioning has failed.

- ☞ The downstream pipe system is leaking.
- ⇒ Check for tightness, see section "Tightness Test", page 9.

### Error Tightness test

The tightness test has failed.

- ☞ The downstream pipe system is leaking.
- ⇒ Check for tightness, see section "Tightness Test", page 9.

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#### Error Pressure sensor

Faulty pressure sensor signal.

- ☞ Cables have been reversed on terminals 3 and 2 (+/- on pressure sensor).
- ⇒ *Check wiring plan and correct, see section “Wiring”, page 5.*
- ☞ Cables defective (lead fracture).
- ⇒ *Check cables.*
- ☞ Faulty pressure sensor signal.
- ⇒ *Open laboratory control unit LCU. The measuring instrument is to be connected in series between the “+” or “-” terminal and the disconnected cable. Measure the current signal: correct current signal range 4–20 mA (4 mA = 0 mbar / 20 mA = 100 mbar)*
- ☞ Laboratory control unit LCU defective.
- ⇒ *Remove the unit and return it to the manufacturer for inspection.*

#### Error Gas deficiency

Inlet pressure  $p_{e \min}$  of the selected inlet pressure range has not been reached.

- ☞ The gas supply is shut off.
- ⇒ *Open the ball valve upstream of the VCL or the upstream valve.*
- ☞ The inlet pressure has fallen below the inlet pressure  $p_{e \min}$  during operation.
- ⇒ *Correct inlet pressure  $p_e$ .*
- ⇒ *Replace the LPG canister if there is no longer enough LPG available.*

#### Error Excess pressure

Inlet pressure  $p_{e \max}$  of the selected inlet pressure range has been exceeded.

- ☞ The inlet pressure has exceeded the inlet pressure  $p_{e \max}$  during operation.
- ⇒ *Correct inlet pressure  $p_e$  using on-site pressure regulator.*

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#### System fault no. 3

CRC checksum error has been detected in the program code.

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
#### System fault no. 5

Error in monitoring process detected.

- ⇒ *Remove the unit and return it to the manufacturer for inspection.*

#### System fault no. 6

Program run error (step error) detected.

- ⇒ *Remove the unit and return it to the manufacturer for inspection.*

#### System fault no. 7

Timing error detected (step time exceeded).

- ☞ Check the downstream volume. This must be between 0.2 and 60 l.
- ⇒ *Please contact our hotline, see section “Product Liability and Warranty”, page 2.*

#### System fault no. 8

Error in communication with monitoring processor detected.

- ⇒ *Remove the unit and return it to the manufacturer for inspection.*

#### System fault no. 9

Relay error detected.

- ☞ Internal system fault.
- ⇒ *Remove the unit and return it to the manufacturer for inspection.*

#### Error Emergency Stop button pressed

Emergency Stop button not connected:

- ☞ The Emergency Stop button is not connected and bridge between terminals 1 and 2 is missing.
  - ☞ Replace bridge between terminals 1 and 2.
- Emergency Stop button connected:
- ☞ Unlock Emergency Stop button

#### Error Display does not light up

The LCU does not respond to “Voltage on” and “Key-operated switch on”, even though voltage is applied.

- ☞ Disconnect the system from the power supply, open LCU and check fuse F1 (T 3.15 AH, 250 V, 5 x 20 mm).
- ☞ Remove the unit and return it to the manufacturer for inspection.

### 10 Accessories

#### Pressure test fitting

To vent the pipe system upstream of the laboratory safety valve VCL in a controlled way, a pressure test fitting can be installed at the inlet of the first valve V1 of the VCL. The inlet pressure  $p_e$  can be displayed with a pressure gauge connected to this pressure test fitting. To enable external tightness checks, nitrogen is supplied via the pressure test fitting.

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**Enclosure:** IP 54

**Housing color:** RAL 7035 light gray

**Configurable, potential-free signaling relay:**

max. 0.1 A, 230 VAC.

3 valve outputs for VCL

Electric current: 315 mA

1 valve output for VCC

Electric current: 500 mA

**Voltage:** 115 VAC, 230 VAC

LCD display for status and faults: 2 x 16 characters.

Activation lock by key-operated switch  
(except LCU..M / LCU..S).

### VCL/VCC

#### Types of gas

Natural gas, LPG (gaseous) and air. The gas must be dry in all temperature conditions and must not contain condensate.

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Valve seal: NBR

Connection flange with internal thread:

Rp acc. to ISO 7-1

### VCL pressure sensor

Prefabricated electrical connection:

plug with socket acc. to EN 175301-803-C

**Socket:** GDSN 307 black

**Enclosure:** IP 65

**Number of pins:** 2 + shield

**Cable gland:** PG 7

**Cable type:** length 5 m, LIYCY, max. 2 x 0.75 mm<sup>2</sup>,  
shielded

Cable ends and shield with wire end ferrules prepared for connection to the LCU. Connect the shield to the LCU (not to the plug).

### Service life

The Pressure Equipment Directive (PED) and the Energy Performance of Buildings Directive (EPBD) demand regular checks on and maintenance of heating systems in order to ensure a high degree of efficiency, smooth operation and functional reliability in the long term.

The design life of the product, hereinafter referred to simply as the “service life”, is based on the relevant standards. You can find further explanations in the applicable rules and regulations and on the afecor website ([www.afecor.org](http://www.afecor.org)).

This service life information presupposes that the product is used in accordance with these operating instructions.

The product must be serviced at regular intervals. Once the end of the service life is reached, all safety-related functions must be checked as specified under “Maintenance”, page 14.

If the product passes the specified function tests, you can continue to use it until the next scheduled maintenance operation. At this point, these tests must be repeated.

If the product fails one of the specified tests, it must be replaced immediately.

This procedure applies to heating systems. For thermo-processing equipment, observe national regulations.

Service life for VAS (from date of manufacture) in accordance with EN 13611, EN 161:

Continuous operation in the upper ambient temperature range accelerates the aging of the elastomer materials and reduces the service life (please contact manufacturer).